Package 'rethomics'

June 4, 2015

2 ethogramPlot

sieepAnnotation	 ٠	•	 ٠	٠	•	 	٠	٠	٠	 •	٠	٠	٠	•	 ٠	٠	•	٠	•	•	•	•	•	П	
virtualBeamCrossClassif						 																		12	9

Index 14

days

Trivially converts days to seconds

Description

Trivially converts days to seconds

Usage

days(x)

Arguments

Χ

number of seconds

Value

the corresponding number of seconds

See Also

hours mins

ethogramPlot

Displays the temporal and inter-individual average of a variable of interest.

Description

This function produces a graph where y is the variable of interest and x, the time. It can be used to visualise temporal trends per groups of conditions. The response variable, y, is grouped by time windows of defined size.

```
ethogramPlot(y, data, condition = NULL, summary_time_window = mins(30),
    normalise_var_per_id = FALSE, error_bar = NULL)
```

fetchPsvResultFiles 3

Arguments

Value

A ggplot object that can be plotted directly, or modified

See Also

oveviewPlot To show per-individual patterns

Examples

```
# Load sample data
data(sleep_validation)
my_data <- sleep_validation[,sleepAnnotation(.SD),by=key(sleep_validation)]
# No condition
p <- ethogamPlot(activity,my_data)
print(p)
# We make a dummy condition, males are in every even region
my_data[,sex:=ifelse(region_id %% 2==0,"m","f")]
p <- ethogamPlot(activity,my_data,condition = sex)
print(p)
p <- ethogamPlot(activity,my_data,condition = sex,error_bar="sm")
# p is simply a ggplot object, so we can change things:
print(p + labs(title="MY own title"))</pre>
```

fetchPsvResultFiles

Query files from a PSV data directory according to the date of the experiment and the device which acquired the data.

Description

This function is designed to list and select experiemental files. In general, end-users will want to retreive path to their experimental files according to the date and ID of the video monitor without having to undertand the underlying directory structure.

```
fetchPsvResultFiles(result_dir, query = NULL)
```

4 fetchPsvResultFiles

Arguments

```
result_dir The location of the result directory (i.e. the folder containing all data).

query An optionnal query formated as a dataframe (see details).
```

Details

The optionnal argument query is expected to be a table where every row maps an experiment. In many respects, it is similar to the what argument in loadPsvData. The only difference is that it does not have a path column. Instead, it must contain two columns:

- time The date and time when the experiment started formated either as 'yyyy-mm-dd' or 'yyyy-mm-dd_hh:mm:ss'. In the former case, there may be several matching experiments to a single time (starting the same day). When this happends, *only the last* is returned, and a warning message is displayed.
- machine_name The name of the machine that acquired the data.

The result is meant to be used directly, as the what argument, by loadPsvData (see examples).

Value

The query extended with the requested paths. When query is not specified, the function returns a table with all available files.

Note

PSV stores the data in a hard-coded directory structure/root_dir/machine_id/machine_name/datetime/file.db:

- machine_id In principle, a universally unique identifier of the acquisition device.
- machine_name A human friendly name for acquisition device. In practice, this is expected to be unique within laboratory.
- datetime The date (and optionnally the time) of the begining of an experiment

Examples

```
## Not run:
# This is where I store the data on my computer
MY_DATA_DIR <- "/data/psv_results/"

query <- data.table(date="2015-06-02", machine_name=c("GGSM-001","GGSM-003"),region_id = rep(1:10,each=2))
print(query)
map <- fetchPsvResultFiles(MY_DATA_DIR, query)
dt <- loadPsvData(map)

## End(Not run)</pre>
```

hours 5

hours

Trivially converts hours to seconds

Description

Trivially converts hours to seconds

Usage

hours(x)

Arguments

Х

number of seconds

Value

the corresponding number of seconds

See Also

days mins

loadDAMFile

Read a text file formatted as DAM2 into a single data table.

Description

This function is used to load data from DAM2 devices as a data.table.

Usage

```
loadDAMFile(FILE, channels = NULL, min_time = 0, max_time = Inf,
interval = 60)
```

Arguments

FILE the name of the input file.

min_time exclude data before min_time (in seconds). This time is relative to the start of

the experiement.

max_time exclude data after max_time (in seconds). This time is relative to the start of the

experiement.

Value

If rois has only one element, a dataframe. Otherwise, a list of dataframes (one per ROI)

6 loadMetaData

Note

Analysis of many long (sevaral days) recording can use a large amount of RAM. Therefore, it can sometimes be advantageaous to load an process ROIs one by one.

See Also

loadMetaData To display global informations about the experiment.

Examples

```
## Not run:
FILE <- "Monitor53.txt"
out <- loadDAMFile(FILE)
#histogram of x marginal distribution
hist(out[roi_id == 1, x], nclass=100)

## End(Not run)
## Not run:
# More realistec example where we have experiemental conditions, and
we want to resample data at 1.0Hz.
# First, the conditions:
conditions <- cbind(roi_id=1:32, expand.grid(treatment=c(T,F), genotype=LETTERS[1:4]))
print(conditions)

## End(Not run)</pre>
```

loadMetaData

Retreive metadata from a result file.

Description

This function is used to obtain metadata – such as 'time and date of the experiment', 'acquisition device', 'version of the software' and such– embeded in a result file generated by PSV.

Usage

```
loadMetaData(FILE)
```

Arguments

FILE

the name of the input file.

Value

A list containing fields for metadata entries

See Also

loadPsvData to load raw data.

IoadPsvData 7

Examples

```
## Not run:
FILE <- "result.db"
out <- loadMetaData(FILE)
names(out)
## End(Not run)</pre>
```

loadPsvData

Read data from a result file.

Description

This function is used to convert all the information contained in a result file generated by PSV (i.e a .db file) into an R 'data.table'.

Usage

```
loadPsvData(what, min_time = 0, max_time = Inf, reference_hour = NULL,
  verbose = TRUE, FUN = NULL, ...)
```

Arguments

what	an object describing which $file(s)$ to load and, optionally, associated variables/conditions (see details).
min_time	exclude data before min_time (in seconds). This time is relative to the start of the experiement.
max_time	exclude data after max_time (in seconds). Also relative to the start of the experiement.
reference_hour	the hour, in the day, to use as t_0 reference. When unspecified, time in the output is relative to the start of the experiment.
verbose	whether to print progress (a logical).
FUN	an optionnal function to transform the data from each 'region' (i.e. a data.table) immedidatly after is has been loaded.
	extra arguments to be passed to FUN

Details

what can be one of two objects:

- A character vector. In which case, it is assumed that each element is the path to a different file to load.
- A dataframe. The dataframe *must* have a column named 'path'. The path basename will be used as a unique identifier for a specific experiement (experiment_id). Arbitrary column can be added to map experimental condition to file name. In addition, the dataframe can have a column named region_id. When defined, only the specified combinations of path and region_id will be loaded. This allows to map additionnal conditions (data frame columns) to specific regions/files. When additionnal conditions are provided, they will result in creation of custom columns in the output of this function.

8 loadPsvData

Value

A data.table where every row is an individual measurment. That is a position at a unique time (t) in a unique region (region_id), and from a unique result file/experiment (experiment_id). The time is expressed in seconds. Distance units (e.g. xy position, height/width) are expressed as a fraction of the with of the region they originate from.

See Also

loadMetaData To display global informations about a specific file.

Examples

```
# First of all, let us load files from the data sample with the package.
# Most likely, you will already have your own data files.
sample_files <- c("validation.db", "sample_1.db", "sample_2.db")</pre>
# Extract the files in your computer
paths <- sapply(sample_files, loadSampleData)</pre>
# Now, `paths` is just a vector of file names:
print(paths)
##################
##################
# Case 1: load ALL REGIONS from a SINGLE FILE
validation_data_file <- paths[1]</pre>
# `validation_data_file` is simply the path to the db file in your computer
dt <- loadPsvData(validation_data_file)</pre>
print(dt)
################
# Case 2: load ALL REGIONS from MULTIPLE FILES
# we pass all the files we want to load as thw `what` argument
dt <- loadPsvData(paths)</pre>
# Note the column `experiment_id` in dt. It tells us which file/experiement
# each measurment originates from.
print(dt)
################
# Case 3: load ALL ROIS from MULTIPLE FILES AND add CONDITIONS
# Let us imagine that each file/experiement
# was acquired under different experiemental condition.
# We can encode this information in a 'master-table' (data frame) in which a column
# named `path` maps experiemental condition(s).
#For instance 3 different treatments:
master_table <- data.frame(path=paths, treatment=c("control", "drug_A", "drug_B"))</pre>
# Let us check our table:
print(master_table)
# The table looks OK, we load the actual data
dt <- loadPsvData(master_table)</pre>
# Note that `dt` now contains a column for your treatment.
print(colnames(dt))
# This makes it easier to perform things such as average per treatment.
print(dt[,list(mean_x = mean(x)),by="treatment"])
################
# Case 4: load SELECTED REGIONS from MULTIPLE FILE, WITH CONDITIONS
# Sometimes, different regions contain for different conditions.
# If the master table has a column named `region_id`, only the specified regions will be returned.
# Let us assume that we want to replicate case 3, but, now, we load only the first 20 regions.
master_table <- data.frame(path=paths, treatment=c("control", "drug_A", "drug_B"), region_id=rep(1:20,each=</pre>
```

IoadSampleData 9

loadSampleData

Retreive sample/example data contained in this package.

Description

This function is intended to be used in order to try and test this package with a reproducible set of raw db data.

Usage

```
loadSampleData(names = NULL)
```

Arguments

names

The name of the samples to be loaded. When names is NULL, the function returns the list of available samples.

See Also

loadPsvData to obtain raw experimental data.

 $\begin{tabular}{ll} max Velocity Classifier & Define whether an animal is moving according to its subpixel velocity. \\ & It requires a variable named $xy_dist_log10x1000$ in the .db file. \\ \end{tabular}$

Description

Define whether an animal is moving according to its subpixel velocity. It requieres a variable named xy_dist_log10x1000 in the .db file.

```
maxVelocityClassifier(data, velocity_threshold = 0.005)
```

10 overviewPlot

Arguments

data the data.table containing behavioural features used for movement classification. velocity_threshold

See Also

sleepAnnotation to apply this function to all subsequent time windows.

mins

Trivially converts minutes to seconds

Description

Trivially converts minutes to seconds

Usage

mins(x)

Arguments

х

number of seconds

Value

the corresponding number of seconds

See Also

days hours

overviewPlot

Displays, per individual, the temporal average of a variable of interest.

Description

This function produces a tiled representation in which every row represents one individual (i.e. from a unique combination of region and experiement). The x axis represents time in days. The values of the variable of interest are represented by different colour intensity.

```
overviewPlot(y, data, condition = NULL, summary_time_window = mins(30),
    normalise_var_per_id = TRUE)
```

sleepAnnotation 11

Arguments

Value

A ggplot object that can be plotted directly or modified.

See Also

ethogramPlot To show trend by aggregating individuals over time.

Examples

```
# Load sample data
data(sleep_validation)
my_data <- sleep_validation[,sleepAnnotation(.SD),by=key(sleep_validation)]
# No condition
p <- overviewPlot(activity,my_data)
print(p)
# We make a dummy condition, males are in every even region
my_data[,sex:=ifelse(region_id %% 2==0,"m","f")]
p <- overviewPlot(activity,my_data,condition = sex)
print(p)
# p is simply a ggplot object, so we can change things:
print(p + labs(title="MY own title"))</pre>
```

sleepAnnotation

Determine whether an animal is asleep

Description

This function uses a motion classifier to first decide whether an animal is moving during a given time window. Then, it defines sleep as contiguous immobility for a minimal duration.

```
sleepAnnotation(data, time_window_length = 10, min_time_immobile = 60 * 5,
  motion_classifier_FUN = maxVelocityClassifier, ...)
```

12 virtualBeamCrossClassif

Arguments

```
the data (i.e a data.table) from a single ROI. It must contain, at least, the columns 't', 'x' and 'y'.

time_window_length

The number of seconds to be used by the motion classifier. This corresponds to the period of the output data.

min_time_immobile

the minimal duration (in s) after which an immobile animal is scored as 'asleep'.

motion_classifier_FUN

the function used to classify movement.

... extra arguments to be passed to motion_classifier_FUN
```

Value

A data table similar to data with additionnal variables/annotations (i.e. 'moving', 'asleep').

Note

The resulting data will only have one data point every time_window_length seconds.

See Also

loadPsvData to load data and optionnaly apply analysis on the fly.

Examples

```
# We load samples from the package data
file <- loadSampleData("validation.db")
# We would like only ROI #2 from this file
map <- data.frame(path=file, roi_id=2)
dt <- loadPsvData(map)
sleep_dt <- sleepAnnotation(dt)
# A more likely scenario: we load ROIs 5 to 10 and
# apply sleep analysis in combination with loadPsvData.
# this means we apply the function to all rois immediatly after they are being loaded.
map <- data.frame(path=file, roi_id=5:10)
dt <- loadPsvData(map,FUN=sleepAnnotation)</pre>
```

virtual Beam Cross Classif

Define whether an animal is moving. This is achieved by computing the number of crossed of a "virtual beam" in the middle of the ROI (i.e. at x=0.5). This emulate the type of data generated by DAM2.

Description

Define whether an animal is moving. This is achieved by computing the number of crossed of a "virtual beam" in the middle of the ROI (i.e. at x=0.5). This emulate the type of data generated by DAM2.

virtualBeamCrossClassif 13

Usage

virtual Beam Cross Classif (data)

Arguments

data

the data.table containing behavioural features used for movement classification.

See Also

 ${\tt maxVelocityClassifier}$ to defince movement by maximum velocity, which is more accurate, instead.

Index

```
days, 2, 5, 10
ethogramPlot, 2, 11
fetchPsvResultFiles, 3
hours, 2, 5, 10
loadDAMFile, 5
loadMetaData, 6, 6, 8
loadPsvData, 4, 6, 7, 9, 12
loadSampleData, 9
maxVelocityClassifier, 9, 13
mins, 2, 5, 10
overviewPlot, 10
overviewPlot, 3
sleepAnnotation, 10, 11
virtualBeamCrossClassif, 12
```